

Product Datasheet

Neuregulin 3 Immunizing Peptide (orb17381)

Description	NRG3 Immunizing peptide
Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	This peptide is used as immunogen to produce Neuregulin 3 antibody (orb19535)
Protein Sequence	RNEIQRDSALTK
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)