

Product Datasheet

CHFR antibody [C3], C-term (orb1741723)

Catalog Number	orb1741723
Category	Antibodies
Description	CHFR antibody [C3], C-term
Target	checkpoint with forkhead and ring finger domains
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Buffer/Preservatives	0.01% Thimerosal. 1XPBS (pH7), 40% Glycerol
Purification	Purified by antigen-affinity chromatography.
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human CHFR. The exact sequence is proprietary.
UniProt ID	Q96EP1
MW	73
Tested applications	WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution range	WB: 1:500-1:3000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.
Application notes	Positive control: H1299 , NT2D1
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)