

Product Datasheet

Human ARTN Protein, hFc Tag (orb1743264)


Catalog Number	orb1743264
Description	Recombinant Human ARTN Protein with C-terminal human Fc tag
Reactivity	Human
Target	ARTN
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-Human Fc Tag
Note	For research use only
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Source	E.coli
MW	The protein has a predicted molecular mass of 38.1 kDa after removal of the signal peptide. The apparent molecular mass of ARTN-hFc is approximately 35-55 kDa due to glycosylation.
Uniprot ID	Q5T4W7
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)