

Product Datasheet

APOBEC3D Immunizing Peptide (orb17774)

Catalog Number	orb17774
Category	Proteins
Description	APOBEC3D Immunizing peptide
Buffer/Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Protein Sequence	ATDSQETRPGRA
Application notes	This peptide is used as immunogen to produce APOBEC3D antibody (orb13204)
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)