

Product Datasheet

Anti-HA-tag Antibody-FITC labeled (orb1786244)

Catalog Number	orb1786244
Category	Antibodies
Description	Rabbit polyclonal antibody to HA-tag-FITC labeled
Species/Host	Rabbit
Conjugation	Unconjugated
Concentration	1 mg/ml
Buffer/Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of HA-tag.
Tested applications	FC, IF, WB
Dilution range	WB (1/2000 - 1/5000), IF (1/500 - 1/2000), FC (1/1000 - 1/4000)
Antibody Type	Secondary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)