

Product Datasheet

Rabbit Anti-Goat IgG (H&L) - Cy5.5 (orb1786338)

Catalog Number	orb1786338
Category	Antibodies
Description	Rabbit Polyclonal Secondary Antibody to Goat IgG (H&L) Cy5.5 labeled
Species/Host	Rabbit
Conjugation	Cy5.5
Reactivity	Goat
Concentration	1 mg/ml
Buffer/Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Goat IgG
Tested applications	ELISA, FC, IF, IH
Dilution range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)