

Product Datasheet

Goat Anti-Mouse IgM mu chain - FITC (orb1786684)

Catalog Number	orb1786684
Category	Antibodies
Description	Goat Polyclonal Secondary Antibody to Mouse IgM mu chain FITC labeled
Species/Host	Goat
Conjugation	FITC
Reactivity	Mouse
Concentration	1 mg/ml
Buffer/Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Mouse IgM mu chain
Tested applications	FC, IF, IH
Dilution range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000)
Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)