

Product Datasheet

Rabbit Anti-Rat IgG (H&L) - FITC (orb1786693)

Catalog Number	orb1786693
Category	Antibodies
Description	Rabbit Polyclonal Secondary Antibody to Rat IgG (H&L) FITC
Species/Host	Rabbit
Conjugation	FITC
Reactivity	Rat
Concentration	1 mg/ml
Buffer/Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Rat IgG
Tested applications	FC, IF
Dilution range	IF (1/250 - 1/2000), FC (1/250 - 1/2000)
Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)