

Product Datasheet

Goat Anti-Rabbit IgG (H&L) - Biotin (orb1786767)

| | |
|-----------------------------|--|
| Catalog Number | orb1786767 |
| Category | Antibodies |
| Description | Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) Biotin |
| Species/Host | Goat |
| Conjugation | Biotin |
| Reactivity | Rabbit |
| Concentration | 1 mg/ml |
| Buffer/Preservatives | 0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Immunogen | Rabbit IgG |
| Tested applications | ELISA, IH, WB |
| Dilution range | E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000) |
| Antibody Type | Secondary Antibody |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)