

Product Datasheet

Rabbit Anti-Goat IgG (H+L) Unconjugated (orb1786896)

Catalog Number	orb1786896
Category	Antibodies
Description	Affinity purification
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Goat
Concentration	1 mg/ml
Buffer/Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Goat IgG (H+L)
Tested applications	ELISA, IHC, WB
Dilution range	EIA and Western blots: 1:5,000-1:100,000. Immunohistochemistry: 1:500-1:5,000. Optimal working dilutions must be determined by end user.
Antibody Type	Secondary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)