

Product Datasheet

PGP9.5 / Uchl1 (pan-Neuronal Marker) Antibody - With BSA and Azide (orb1790130)

| | |
|----------------------------|--|
| Catalog Number | orb1790130 |
| Category | Antibodies |
| Description | Mouse Monoclonal Antibody [Clone SPM575] |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Conjugation | Unconjugated |
| Reactivity | Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Sheep, Zebrafish |
| UniProt ID | P09936 |
| Tested applications | IF, WB |
| Antibody Type | Primary Antibody |
| Clone Number | SPM575 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)