

Product Datasheet

CD56 antibody (C-term) (orb1790244)

Catalog Number	orb1790244
Category	Antibodies
Description	Purified Rabbit Polyclonal Antibody (Pab)
Target	NCAM1 (HGNC:7656)
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Rabbit Ig
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Immunogen	This antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between amino acids from human.
UniProt ID	P13591
MW	94574 Da
Tested applications	WB
Dilution range	WB - 1:1000
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note


For research use only

Expiration Date

12 months from date of receipt.



All lanes: Anti-CD56 antibody (C-term) at 1:1000 dilution. Lane 1: mouse brain lysate. Lane 2: rat brain lysate. Lane 3: U-87 MG whole cell lysate. Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated E42 at 1/15000 dilution. Observed band size: 120-220KDa. Blocking/Dilution buffer: 5% NFDM/TBST.



All lanes: Anti-CD56 antibody (C-term) at 1:1000 dilution. Lane 1: mouse brain lysate. Lane 2: rat brain whole cell lysate. Lane 3: U-87 MG whole cell lysate. Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated E42 at 1/15000 dilution. Observed band size: 150KDa. Blocking/Dilution buffer: 5% NFDM/TBST.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558