

Product Datasheet

Anti-Hu CD85j Purified Low Endotoxin (orb179908)

Description	Mouse monoclonal antibody to CD85j.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IP, WB
Immunogen	Hairy cell leukaemia cells
Target	CD85j
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Clone Number	GHI/75
Antibody Type	Primary Antibody
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q8NHL6

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Entrez 10859


Expiration Date 12 months from date of receipt.



Flow cytometry multicolor surface staining of human lymphocytes stained using anti-human CD85j (GHI/75) purified antibody (low endotoxin, concentration in sample 1 µg/ml) GAM APC and anti-human CD19 (LT19) PE antibody (20 µl reagent / 100 µl of peripheral whole blood).



Separation of human CD85j positive B cells (red-filled) from neutrophil granulocytes (black-dashed) in flow cytometry analysis (surface staining) of human peripheral whole blood stained using anti-human CD85j(GHI/75) purified antibody (low endotoxin, concentration in sample 1 µg/ml) GAM APC.



Flow cytometry surface staining pattern of human peripheral blood stained using anti-human CD85j(GHI/75) purified antibody (low endotoxin, concentration in sample 1 µg/ml) GAM APC.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)