

Product Datasheet

Chloramphenicol antibody (orb1804811)

Catalog Number	orb1804811
Description	Chloramphenicol antibody
Species/Host	Rabbit
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	chloramphenicol-HMS-BSA
Target	Chloramphenicol
Preservatives	10mM sodium azide. PBS, pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
MW	1kDA
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)