

Product Datasheet

Anti-ERBB1 / EGFR / HER1 Reference Antibody (tomuzotuximab) (orb1805550)

Catalog Number	orb1805550
Category	Antibodies
Description	Anti-ERBB1 / EGFR / HER1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light chain type is hukappa. It has a predicted MW of 145.5 kDa.
Target	ERBB1 / EGFR / HER1
Clonality	Monoclonal
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	100 mM Proline/Acetic acid 20 mM Arginine pH 5.0
Purity	>95%
MW	145.5 kDa
Tested applications	ELISA, FA, FACS, Kinetics
Endotoxins	1 EU/mg
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Anti-ERBB1 / EGFR / HER1 Reference Antibody (tomuzotuximab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%.



The purity of Anti-ERBB1 / EGFR / HER1 Reference Antibody (tomuzotuximab) is more than 95% ;determined by SEC-HPLC.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)