

Product Datasheet

Rat RANK ELISA Kit (orb180647)

Catalog Number	orb180647
Category	Assays and Kits
Description	Rat RANK ELISA Kit
Target	TNF receptor superfamily member 11A
Reactivity	Rat
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q30-P213
UniProt ID	F1M8Z6
MW	67607 Da
Tested applications	ELISA
Specificity	Natural and recombinant rat RANK
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Standard	10 ng/tube
Range	62.5 pg/ml - 4,000 pg/ml
Sensitivity	10 pg/ml
Sample Types	cell culture supernatants and serum.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Rat RANK PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)