

Product Datasheet

Recombinant human ADAM17 protein, His (HEK293) (orb1808160)

Catalog Number	orb1808160
Description	Recombinant human ADAM17 protein, His (HEK293)
Target	Disintegrin and metalloproteinase domain-containing protein 17
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	358-649/824
Purity	> 90% as determined by SDS-PAGE
MW	34 kDa
Expiration Date	6 months from date of receipt.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)