

## Product Datasheet

### Rat IGFBP2 ELISA Kit (orb1809056)

<b>Catalog Number</b>	orb1809056
<b>Category</b>	Assays and Kits
<b>Description</b>	The Rat IGFBP2 ELISA Kit features a dynamic range of 3.12-200ng/mL and is suitable for Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. This product offers a sensitivity of 1.34ng/mL.
<b>Reactivity</b>	Rat
<b>Range</b>	3.12-200ng/mL
<b>Sensitivity</b>	1.34ng/mL
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)