

Product Datasheet

Human ACE2 ELISA Kit (orb1809378)

Catalog Number	orb1809378
Category	Assays and Kits
Description	The Human ACE2 ELISA Kit features a dynamic range of 15.6-1000pg/mL and is suitable for Serum, plasma or other biological fluids. This product offers a sensitivity of 7.7pg/mL.
Reactivity	Human
Range	15.6-1000pg/mL
Sensitivity	7.7pg/mL
Sample Types	Serum, plasma or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)