

## Product Datasheet

# Vitronectin Recombinant Rabbit Monoclonal Antibody (orb1816904)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Vitronectin Recombinant Rabbit Monoclonal Antibody                                                                                      |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | IF, IHC-Fr, IHC-P, WB                                                                                                                   |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human Vitronectin                                                                         |
| <b>Target</b>              | VTN                                                                                                                                     |
| <b>Preservatives</b>       | 1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.02% Proclin300.                                                            |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                  |
| <b>Concentration</b>       | 1mg/ml                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Recombinant                                                                                                                             |
| <b>Clone Number</b>        | 7A1                                                                                                                                     |
| <b>Antibody Type</b>       | Recombinant Antibody                                                                                                                    |
| <b>MW</b>                  | 52 kDa                                                                                                                                  |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

|                 |                                                              |
|-----------------|--------------------------------------------------------------|
| Uniprot ID      | P04004                                                       |
| Dilution Range  | WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:400-800, IF=1:50-200 |
| Expiration Date | 12 months from date of receipt.                              |

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)