

Product Datasheet

F(ab')₂ Fragment Goat Anti-rabbit IgG H&L, FITC conjugated (orb1817038)

Catalog Number	orb1817038
Category	Antibodies
Description	F(ab') ₂ Fragment Goat Anti-rabbit IgG H&L, FITC conjugated
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	FITC
Reactivity	Rabbit
Buffer/Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Purification	Affinity purified, nonspecific adsorbed
Immunogen	Native rabbit IgG
Tested applications	FC, ICC, IF
Dilution range	IF=1:200-1000, ICC/IF=1:100-1000, Flow-Cyt=1:100-1000
Antibody Type	Secondary Antibody
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)