

Product Datasheet

CTLA-4 / CD152 Reference Antibody (zalifrelimab) (orb1817627)

Catalog Number	orb1817627
Category	Antibodies
Description	Recombinant Antibody
Clonality	Monoclonal
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative. This antibody is purified through a protein A column.
Concentration	1mg/ml
Protein Sequence	CTLA-4 / CD152
UniProt ID	P16410
MW	145 Da
Tested applications	FC, Kinetics
Antibody Type	Primary Antibody
Expression System	CHO Cell
Endotoxins	0.001EU/μg, determined by LAL method.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only



The purity of Anti-CTLA-4 / CD152 Reference Antibody (zalifrelimab) is more than 95%, determined by SEC-HPLC.



Anti-CTLA-4 / CD152 Reference Antibody (zalifrelimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558