

Product Datasheet

Anti-IL-1RL1 / ST2 / IL-33R Reference Antibody (melrilimab) (orb1817882)


Description	Recombinant Antibody
Reactivity	Human, Monkey, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Protein Sequence	IL-1RL1 / ST2 / IL-33R
Isotype	IgG2SA
Clonality	Monoclonal
Antibody Type	Primary Antibody
MW	145.2
Uniprot ID	Q01638
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



The purity of Anti-IL-1RL1 / ST2 / IL-33R Reference Antibody (melrilimab) is more than 95%, determined by SEC-HPLC.



Anti-IL-1RL1 / ST2 / IL-33R Reference Antibody (melrilimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)