

Product Datasheet

Apo Active FITC - Caspase 3 Detection Kit (orb1821178)

Catalog Number	orb1821178
Category	Assays and Kits
Description	This kit is highly specific for active human and murine enzyme. The kit utilizes a rabbit affinity purified polyclonal antibody against aa 163-175 of murine caspase 3. Visualization of fixed samples is via a FITC-labeled goat anti-rabbit antibody.
Assay Type	Detection
Storage	2-8°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)