

Product Datasheet

MRSA (PBP2a) Antibody (FITC) (orb1821257)

Catalog Number	orb1821257
Category	Antibodies
Description	MRSA (PBP2a) Antibody (FITC)
Target	MRSA (PBP2a)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a
Conjugation	Unconjugated
Reactivity	Bacteria
Concentration	1mg/ml
Buffer/Preservatives	Purified antibody conjugated to Fluorescein using Innova Biosciences Lightning-Link , supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Purification	Protein A purified
Immunogen	Native Staphylococcus aureus (methicillin resistant).
Tested applications	ICC
Dilution range	As detection antibody in ELISA we suggest a range of 1/5,000 – 1/25,000. Optimal antibody dilution should be determined by user laboratory since it may be dependent on the detection system used in both ELISA and IF
Application notes	SP2/0

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	AC10
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)