

Product Datasheet

Carcino Embryonic Antigen CEA Monoclonal Antibody [CB30] - 20ul-Trial Size (orb1821402)

Description	Carcino Embryonic Antigen CEA Monoclonal Antibody [CB30]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Full length native protein (purified) (Human)
Target	Carcino Embryonic Antigen CEA
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	CB30
Uniprot ID	P06731
Entrez	1048
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at WB 1–µg/ml MCF7, 0.003–µg/ml human colorectal carcinoma, ICC/IF 1:50-1:100, Flow Cyt 1:50-1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)