

## Product Datasheet

### Horse Genomic DNA (orb1821417)

<b>Catalog Number</b>	orb1821417
<b>Category</b>	Tools
<b>Description</b>	Horse genomic DNA is recommended for genomic analysis e.g. as template in positive control PCRs like for the MeatDetect qPCR Kit Pork+Horse (Halal). The DNA was purified from muscle meat of horse. The DNA fragments have a size of up to 50 kb.
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	200 ng/μl
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles. aliquoting of DNA samples is recommended
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

---

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)