

## Product Datasheet

### GNMB Biosimilar Antibody (orb1826895)

<b>Catalog Number</b>	orb1826895
<b>Category</b>	Antibodies
<b>Description</b>	Anti-GNMB(glembatumumab biosimilar) mAb
<b>Target</b>	GNMB
<b>Clonality</b>	Monoclonal
<b>Species/Host</b>	Human
<b>Isotype</b>	Human IgG2(S378A) - Kappa
<b>Conjugation</b>	Unconjugated
<b>Reactivity</b>	Human
<b>Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
<b>Purification</b>	Purified from cell culture supernatant by affinity chromatography
<b>UniProt ID</b>	<b>Q14956</b>
<b>Tested applications</b>	ELISA, FC
<b>Dilution range</b>	ELISA 1:5000-10000, Flow Cyt 1:100
<b>Antibody Type</b>	Biosimilar Antibody
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**


68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


**Note** For research use only

**Expiration Date** 12 months from date of receipt.

**Anti-GPNMB(glembatumumab biosimilar) mAb ELISA**  
0.2 µg of Human GPNMB, His tagged protein per well



ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human GPNMB Protein, His Tag (orb1743304) can bind Anti-GPNMB (glembatumumab biosimilar) mAb (orb1826895) in a linear range of 3.20–80 ng/ml.



Flow cytometry analysis with 1 µg/mL Anti-GPNMB (glembatumumab biosimilar) mAb (orb1826895) on Expi293 cells transfected with Human GPNMB protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558