

Product Datasheet

LSM4 siRNA (Human) (orb1859750)

Catalog Number	orb1859750
Category	Tools
Description	siRNA to inhibit LSM4 expression using RNA interference
Reactivity	Human
Form/Appearance	Lyophilized powder
Purity	> 97%
UniProt ID	Q9Y4Z0
Tested applications	RNA-I
Application notes	We recommends transfection with 10 nM - 100 nM siRNA 48 to 72 hours prior to cell lysis. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of the tube. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. For example, resuspend one tube of 5 nmol siRNA oligo in 250 µl of DEPC water to get a final concentration of 20 µM.
Source	Synthetic
Storage	Shipped at 4 °C. Store at -20 °C for one year.
Note	For research use only
Entrez	25804

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)