

Product Datasheet

UBE2G2 siRNA (Human) (orb1863248)

Catalog Number	orb1863248
Description	siRNA to inhibit UBE2G2 expression using RNA interference
Reactivity	Human
Tested Applications	RNA-I
Form/Appearance	Lyophilized powder
Storage	Shipped at 4 °C. Store at -20 °C for one year.
Note	For research use only
Application notes	We recommends transfection with 10 nM - 100 nM siRNA 48 to 72 hours prior to cell lysis. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of the tube. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. For example, resuspend one tube of 5 nmol siRNA oligo in 250 µl of DEPC water to get a final concentration of 20 µM.
Purity	> 97%
Source	Synthetic
Uniprot ID	P60604
Entrez	7327
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)