

Product Datasheet

Recombinant Human Myelin P2 protein/PMP2 (orb1906990)

Catalog Number	orb1906990
Category	Proteins
Description	Recombinant Human Myelin P2 protein/PMP2 is produced by E. coli. The target gene encoding S2-V132 is expressed with a His tag at the C-terminus.
Tag	C-His
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
UniProt ID	P02689
Expression System	E. coli
Biological Origin	Human
Expression Region	S2-V132
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)