

Product Datasheet

Recombinant Human IL-15RA & IL-15 fusion protein (Fc) (orb1907226)

Catalog Number	orb1907226
Category	Proteins
Description	Recombinant Human Interleukin-15 Receptor Alpha & Interleukin-15 Fusion Protein is produced by our Mammalian expression system and the target gene encoding Ile31-Asp96&Asn49-Ser162(Asn120Asp) is expressed with a human IgG1 Fc tag at the C-terminus.
Tag	C-hFc
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization.
UniProt ID	P40933, Q13261
Expression System	Human 293 Cells
Biological Origin	Human
Expression Region	Ile31-Asp96 & Asn49-Ser162(Asn120Asp)
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)