

Product Datasheet

Influenza B NP Human Monoclonal Antibody (Detection) (orb1921577)

Catalog Number	orb1921577
Category	Antibodies
Description	Influenza B NP Human Monoclonal Antibody (Detection)
Target	Influenza B
Clonality	Monoclonal
Species/Host	Human
Conjugation	Unconjugated
Buffer/Preservatives	NaN3
Tested applications	EIA
Matched Pairs	Capture: orb1921576; Detection: orb1921577;
Source	HEK Cells
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)