

Product Datasheet

Soluble Growth Stimulation Expressed Gene2 (ST2) (His Tag) Mouse Monoclonal Antibody (Capture) (orb1921753)

Catalog Number	orb1921753
Category	Antibodies
Description	Soluble Growth Stimulation Expressed Gene2 (ST2) (His Tag) Mouse Monoclonal Antibody (Capture)
Clonality	Monoclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Purified,Liquid
Buffer/Preservatives	PBS, pH 7.4
Purity	>95% pure (SDS-PAGE)
Purification	Protein A Chromatography
MW	150 kDa
Tested applications	CLIA
Specificity	Specific for ST2
Matched Pairs	Control:orb1921755; Capture: orb1921753; Detection: orb1921754;

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Source	Mammalian cell
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)