

Product Datasheet

Goat Anti-Rabbit IgG H&L (Cy5) (orb1923984)

Description	Goat Anti-Rabbit IgG H&L (Cy5)
Reactivity	Rabbit
Conjugation	Cy5
Tested Applications	ELISA, FC, ICC, IF, IHC-Fr
Immunogen	Rabbit IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	IF: 1/100-1/1000 FC: 1/500-1/1000 ELISA: 1/5000-1/20000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)