

Product Datasheet

Goat Anti-Mouse IgG H&L (Biotin) (orb1923998)

Catalog Number	orb1923998
Category	Antibodies
Description	Goat Anti-Mouse IgG H&L (Biotin)
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Biotin
Reactivity	Mouse
Form/Appearance	Liquid
Concentration	1 mg/ml
Buffer/Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Purification	Affinity Purified
Immunogen	Mouse IgG
Tested applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Dilution range	WB: 1/5000-1/20000 IF: 1/100-1/1000 ELISA: 1/5000-1/20000
Antibody Type	Secondary Antibody
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)