

Product Datasheet

Loading Control Rabbit pAb Kit Antibody (orb1924029)

Catalog Number	orb1924029
Category	Assays and Kits
Description	Loading Control Rabbit pAb Kit Antibody
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	liquid
Buffer/Preservatives	Beta actin: Rabbit IgG in phosphate buffered saline,pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Beta Tubulin: Rabbit IgG in phosphate buffered saline,pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. GAPDH: Rabbit IgG in phosphate buffered saline,pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography
UniProt ID	P07437, P60709, P04406
MW	Beta actin: 42 kDa; Beta Tubulin: 55 kDa; GAPDH: 36 kDa
Tested applications	ICC, IF, IHC-P, WB
Dilution range	WB: 1/5000-1/10000 IHC: 1/50-1/200 IF: 1/100-1/200
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)