

Product Datasheet

Goat Anti-Rabbit IgG H&L (TRITC) (orb1924911)

Description	Goat Anti-Rabbit IgG H&L (TRITC)
Reactivity	Rabbit
Conjugation	TRITC
Tested Applications	FC, ICC, IF
Immunogen	Rabbit IgG
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	liquid
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	IF: 1/100-1/1000 FC: 1/500-1/1000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)