

Product Datasheet

Anti-Zika NS1 IgG (orb1925042)

Description	This chimeric whole IgG antibody contains murine variable region for Zika NS1 protein binding and human IgG constant region for detection by an anti-human IgG Fc antibody.
Conjugation	Unconjugated
Preservatives	0.1% ProClin300. cell culture medium
Form/Appearance	0.22 um filtered cell culture supernatant
Storage	Upon receipt, store at 2-8°C
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)