

## Product Datasheet

### DCLK1 Antibody (C-term) (orb1927559)

|                      |                                                                                                                                                                              |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number       | orb1927559                                                                                                                                                                   |
| Category             | Antibodies                                                                                                                                                                   |
| Description          | DCLK1 Antibody (C-term)                                                                                                                                                      |
| Clonality            | Polyclonal                                                                                                                                                                   |
| Species/Host         | Rabbit                                                                                                                                                                       |
| Isotype              | Rabbit IgG                                                                                                                                                                   |
| Conjugation          | Unconjugated                                                                                                                                                                 |
| Reactivity           | Human                                                                                                                                                                        |
| Predicted Reactivity | Mouse                                                                                                                                                                        |
| Form/Appearance      | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| UniProt ID           | <b>O15075</b>                                                                                                                                                                |
| MW                   | 82224 Da                                                                                                                                                                     |
| Tested applications  | WB                                                                                                                                                                           |
| Dilution range       | WB - 1:1000                                                                                                                                                                  |
| Specificity          | This DCLK1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 685-712 amino acids of human DCLK1.                                  |
| Antibody Type        | Primary Antibody                                                                                                                                                             |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States


Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

**Note**

For research use only



DCLK1 Antibody (C-term) western blot analysis in MDA-MB231 cell line lysates (35 ug/lane). This demonstrates the DCLK1 antibody detected the DCLK1 protein (arrow).

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)