

Product Datasheet

Boc-GABA-OH (orb1940674)

Catalog Number	orb1940674
Description	Boc-GABA-OH is a PROTAC linker that can be used to synthesize UNC6852, an EED-targeted PROTAC.
Target	VEGFR2
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_9H_{17}NO_4$
Purity	> 98% (HPLC)
MW	203.2
CAS Number	57294-38-9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)