

Product Datasheet

5'-FAM-

GGUGAAAGACGAUGGACAATT,UUGUCCAUCGUCUUUCAC
(orb1941473)

Catalog Number	orb1941473
Category	Tools
Description	This product is a 5'-FAM-labeled double-stranded RNA oligonucleotide with excitation/emission at ~490/520 nm. Such fluorescent probes are widely used in real-time detection for in vitro applications like qPCR assays and for studying RNA interference and cellular localization mechanisms.
Target	Others
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)