

Product Datasheet

CB2R agonist 3 (orb1942001)

Catalog Number	orb1942001
Category	Small Molecules
Description	CB2R agonist 3 is a potent small molecule activator of cannabinoid receptor 2 (CB2R) with an EC50 of 0.37 μ M. Its selectivity for CB2R makes it a valuable research tool for studying the receptor's role in immune modulation, inflammation, and pain in both cellular and animal models.
Target	Cannabinoid Receptor
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)