

Product Datasheet

Erzotabart Antibody (orb1942326)

Catalog Number	orb1942326
Category	Antibodies
Description	Erzotabart is an investigational IgG1-kappa monoclonal antibody targeting human CD38. Its antineoplastic activity, demonstrated in preclinical studies, stems from disrupting CD38-mediated pathways crucial for tumor cell proliferation, metabolism, and immune evasion in oncology research.
Target	CD38
Clonality	Monoclonal
Conjugation	Unconjugated
Concentration	5.18 mg/mL
Purity	95.00% (May vary between batches)
MW	145.01 kDa
Antibody Type	Inhibitor Antibody
Target Areas	CD38
CAS Number	2430792-01-9
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558