

Product Datasheet

Human IgG1 (D265A) kappa, Isotype Control Antibody (orb1942605)

Catalog Number	orb1942605
Description	Human IgG1 (D265A) kappa Isotype Control is a chimeric monoclonal antibody derived from both human and mouse, designed to serve as an isotype control for the human IgG1 (D265A) κ antibody.
Conjugation	Unconjugated
Target	Human IgG1 (D265A) kappa, Isotype Control
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Purity	98.00%
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)