

Product Datasheet

Neuromedin S(rat) TFA (orb1943083)

Catalog Number	orb1943083
Description	Neuromedin S(rat) TFA, a 34-amino acid peptide derivative of rat Neuromedin S, is a neuropeptide sourced from the rat brain that functions as a ligand for the G protein-coupled receptor FM4/TGR-1 [1].
Storage	-20°C
Note	For research use only
Purity	98.00%
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)