

Product Datasheet

Prosaptide Tx14(A) TFA (orb1943378)

Catalog Number	orb1943378
Category	Proteins
Description	Prosaptide Tx14(A) TFA is a prosaposin-derived peptide that potently activates GPR37 and GPR37L1 receptors. It induces ERK1/2 phosphorylation in Schwann cells and is a valuable research tool for studying neuroprotection, myelination, and glioblastoma in both cellular and animal models.
Target	ERK
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)