

Product Datasheet

TGR5 Receptor Agonist 4 (orb1943763)

Catalog Number	orb1943763
Category	Small Molecules
Description	TGR5 Receptor Agonist 4 is a potent and selective small molecule agonist for the bile acid receptor TGR5, demonstrating EC50 values of 2 nM (hTGR5) and 3 nM (mTGR5). It is a valuable research tool for studying TGR5-mediated metabolic functions, including glucose homeostasis and body weight regulation, in both in vitro and in vivo models.
Target	Others, GPCR19
MW	470.95
CAS Number	1872436-42-4
Formula	C ₂₅ H ₂₇ ClN ₂ O ₅
SMILES	<chem>C(OCC1=C(OC)C=CC=C1)(=O)C=2[C@@H](CC(=O)N(C[C@H]3CCCCO3)C2C)C=4C=CC(Cl)=NC4</chem>
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)