

Product Datasheet

TLR4 agonist-1 TEA (orb1943786)

Catalog Number	orb1943786
Category	Small Molecules
Description	TLR4 agonist-1 (TEA, compound 17a) is a potent small-molecule activator of Toll-like Receptor 4, demonstrated to induce MIP-1 β production in RAW 264.7 and MM6 cell lines. This tool compound is used in vitro to study innate immune signaling, inflammation, and TLR4 pathway biology in immunological research.
Target	TLR,Others
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)