

Product Datasheet

4-Methylhistamine (orb1944110)

Catalog Number	orb1944110
Category	Small Molecules
Description	4-Methylhistamine is a potent and selective histamine H4 receptor (H4R) agonist. It is a valuable research tool for investigating the role of H4R in immunology, with applications in in vitro and in vivo studies of inflammation, autoimmune diseases, and cancer.
Target	Histamine Receptor
MW	125.17
CAS Number	36507-31-0
Formula	$C_6H_{11}N_3$
SMILES	<chem>N1=CNC(=C1C)CCN</chem>
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)